

# Haemo-Sol ALT Regular (026-059)



- General, all-purpose cleaner scientifically formulated to perform like Haemo-Sol Regular (026-050).
- Removes protein and organic debris using safer ingredients.
- Maintains efficacy longer than liquid products because surfactants and enzymes don't breakdown or degrade.
- Formulated to be non-corrosive and more environmentally-friendly than hazardous solvents and caustic acids.
- Manufactured under strict cGMP standards and guidelines.

***Haemo-Sol ALT Regular (026-059) is specifically designed to offer the same cleaning benefits of Haemo-Sol Regular (026-050), while being better aligned with the United States Environmental Protection Agency's Safer Chemical Ingredient List. Blended with high-quality chemicals that lift and suspend contaminants in solution, Haemo-Sol ALT Regular leaves items clean, and residue- and contaminant-free after proper rinsing.***

## **Mix Instructions & Directions:**

Dilute ½ ounce (14 grams) of detergent to 1 gallon (3.79 liters) of warm water. Ambient temperature water may be acceptable, depending on type of contaminant being removed. For difficult soils, use hot water. If using an ultrasonic cleaner, use 1/3 to ½ ounce (9 to 14 grams) to 1 gallon (3.79 liters) of water.

**Note:** Hard water full of minerals, such as calcium and magnesium, may prevent the surfactants in Haemo-Sol ALT Regular from working properly. When using hard water, try our Haemo-Sol Enzyme Active or Haemo-Sol Non-Sudsing.

## **Cleaning Method:**

Soak items in Haemo-Sol bath. May use a cloth, brush, sponge, or pad to break-up contaminants. Product can also be used in an ultrasonic, megasonic, or continuous flow rinse tank system. Not intended for use in pressurized sprayers or mechanical cleaners.

## **Tests & Studies:**

Visit [www.haemo-sol.com](http://www.haemo-sol.com) for testing data on and technical information about Haemo-Sol ALT Regular and our other products.

## **Contaminants & Soils Removed:**

Organic by-products, bodily fluids, blood, fats, oils, waxes, tissues, salts, chemicals, pigments, petrolatum, carbomers, alcohols, soil, grit, grime, buffing compounds, slime, grease, particulates, solvents, and silicon oils, to name a few.

## **Used to clean:**

Recommended for glass, stainless steel, silicon wafers, porcelain, ceramic, plastic, rubber, and fiberglass. Used to clean healthcare instruments, laboratory ware, glassware, tissue culture ware, tubing, optical parts, electronic components, stainless steel manufacturing equipment (pharmaceutical, cosmetic, food, beverage, etc.), industrial parts, tanks, electrodes, implantable transmitters, grinders, and more.

## **Industries:**

Pharmaceutical; Laboratory (e.g., Pathology, Hematology, etc.); Veterinary; Dental; Tattoo; Hospital/Medical; Mortuary; Food & Beverage; Agriculture; Meat & Poultry Processing Facilities; and Photovoltaic/Electronics.

